



RECROSS



SCAN CODE FOR INSTRUCTION VIDEOS

MODEL	GUIDE CATHETER COMPATIBILITY	DISTAL SHAFT PROFILE	PROXIMAL SHAFT PROFILE	GUIDEWIRE COMPATIBILITY	USABLE LENGTH
RC144025	5F	2.3F/3.3F	2.6F/3.4F	0.014"	140 cm

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.
 Caution: Federal Law (U.S.) restricts this device to sale by or on the order of a physician.
 In-vitro test data on file at IMDS.
 Bench test data may not be indicative of clinical performance.
 ReCross™ is a trademark of IMDS.


© IMDS 2023 ITM200184 R2.0



Phone: +31(0)50 8200230
 Fax: +31(0)50 8200231

Ceinturbaan Noord 150
 9301 NZ Roden
 The Netherlands
 sales@imds.nl

TOOLS FOR EVERYDAY EFFICIENCY IN COMPLETE REVASCULARISATION



NHANCER Rx
 UNIQUE 6F GUIDING CATHETER COMPATIBILITY SUPPORTING OPTIMAL SIDE BRANCH ACCESS

TRAP
 EFFICIENT, SAFE AND SIMPLE OTW DEVICE EXCHANGE

BLIMP
 DESIGNED TO CROSS

NHANCER PROx
 HIGH PUSH AND RETROGRADE CROSSING PERFORMANCE SUPPORTING YOUR COMPLEX PCI

FLOWGUIDE
 REDUCING PROCEDURAL INDUCED ISCHEMIA

MICRO Rx
 CONVENIENT RAPID EXCHANGE GUIDEWIRE SUPPORT

TAILORED CTO SIDE BRANCH ACCESS

RECROSS

UNIQUE ANTEGRADE PARALLEL WIRING SUPPORT

TWO SIDE PORTS ALLOW WIRE REDIRECTION FROM SUBINTIMAL SPACE



IMDS

INTERVENTIONAL MEDICAL DEVICE SOLUTIONS



UNPARALLELED GUIDEWIRE REDIRECTION OPTIONS





EVERYDAY EFFICIENCY



INTERVENTIONAL MEDICAL DEVICE SOLUTIONS

WWW.IMDS.NL

PRECISE TARGETING OF SIDE BRANCH
OCCLUSION



STRONG FRONTLINE CROSSING AND
PENETRATION WHEN NEEDED


ANTEGRADE
DUAL ACCESS (ADA)

RECCROSS

UNPARALLELED GUIDEWIRE REDIRECTION OPTIONS

360° CIRCUMFERENTIAL ACCESS DUAL LUMEN
ENABLES UNIQUE & EFFECTIVE PROCEDURES

OVAL DISTAL SHAFT
MINIMIZING CROSSING PROFILE




RECCROSS

VERSATILITY FOR EVERYDAY EFFICIENCY


NHANCER Rx



6F GUIDING CATHETER
COMPATIBILITY



SUBINTIMAL GUIDEWIRE
REDIRECTION (SGR)



REMOVABLE STYLET
PROVIDES SIGNIFICANT PUSHABILITY IMPROVEMENT

HEMATOMA
DECOMPRESSION